



Dansac Product Guide

Updated July 2021

**Healthy skin for
a healthy life**

Graham, urostomy since 2014

www.dansac.com.au

Stoma care with a soul

People are more than their stomas. They have real lives. Moving from one day to the next, facing challenges we all face, in addition to the challenges living with a stoma can bring. We know it can't be easy. So it is our job to make sure we do everything we can to make living with a stoma a better experience for everyone.

We start by focusing on the needs of the whole person. Their stoma. Their skin. Their soul. Adjusting to life with a stoma requires a lot of support.

It starts with our products. We design them to work hard to promote a natural skin environment. We know keeping the skin around the stoma healthy goes a long way in enhancing the quality of people's lives.

We make sure our skin-friendly products are easy to use. This not only helps our patients and caregivers, but it also makes our products an easier choice for clinicians.

We help people overcome not only the physical aspects, but also the social and psychological repercussions of getting a stoma. Their lives have undergone a big change and it is important to us that they feel confident in themselves and good in their own skin.

Our dedication to stoma care, from stoma to skin to soul, is the foundation of our brand. It has brought us to where we are today and will lead us to new places in the future.



NovaLife TRE™

NovaLife TRE

NovaLife TRE 1 Closed Flat	10
NovaLife TRE 1 Closed Soft Convex	11
NovaLife TRE 1 Closed Convex	12
NovaLife TRE 1 Open Flat	13
NovaLife TRE 1 Open Soft Convex	14
NovaLife TRE 1 Open Convex	15
NovaLife TRE 2 Flat Wafer	17
NovaLife TRE 2 Convex Wafer	18

NovaLife GX & GX+

NovaLife Closed Pouches

NovaLife GX 1 Closed Flat	21
NovaLife GX 1 Closed Soft Convex	22
NovaLife GX 1 Closed Convex	23

NovaLife Open Pouches

NovaLife GX 1 Open Flat	24
NovaLife GX+ 1 Open Flat	25
NovaLife GX 1 Open Soft Convex	26
NovaLife GX 1 Open Convex	28
NovaLife GX+ 1 Open Convex	29

NovaLife 2 Piece

NovaLife 2 Closed Pouch	30
NovaLife 2 Open Pouch	31
NovaLife GX 2 Flat Wafer	32
NovaLife GX+ 2 Flat Wafer	33
NovaLife GX 2 Convex Wafer	34
NovaLife GX+ 2 Convex Wafer	35

Nova

Nova Closed Pouches

Nova 1 Closed Flat	39
Nova 1 Closed Soft Mouldable	40
Nova 1 Closed Convex	41
Nova 2 Closed Pouch	42

Nova Open Pouches

Nova 1 FoldUp Flat	43
Nova 1 FoldUp X3	44
Nova 1 FoldUp Soft Mouldable	45
Nova 1 FoldUp Convex	46
Nova 1 Infant / Paediatric	47
Nova 2 FoldUp Pouch	48
Nova 2 High Output Pouch	49

Nova Urostomy Pouches

Nova 1 Urostomy X3	50
Nova 1 Urostomy Convex	51
Nova 2 Urostomy Pouch	52

Nova 2 Wafers

Nova 2 Flat Wafer	53
Nova 2 X3 Wafer	54
Nova 2 Soft Mouldable Wafer	55
Nova 2 Convex Wafer	56

Accessories

Dansac TRE Seal, Dansac Seal, Dansac GX-tra Seal	59
Dansac X-tra Strips, Dansac EasiSpray™ Adhesive Remover, Dansac Soft Paste	60
Dansac Skin Lotion Tissue, Dansac Skin Lotion, Dansac Skin Creme	61
Dansac Nodor S, Belts, Clamps	62
Dansac Irrigation Set	63

Index

Index	67
-------------	----



Dansac NovaLife TRE


Shelley, stoma since 2015

Dansac TRE

Technology Explained



- The best skin is **healthy skin**.
- **pH** is a measure of how acidic or alkaline something is. **pH 7** is **neutral**.
- **Healthy skin** has a pH between 4 and 6 which means it is slightly acidic. This is known as the **acid mantle**.
- Maintaining the acid mantle of the skin is essential to help protect its antimicrobial properties and **maintain skin integrity**.
- The acid mantle of the skin **may become weakened** if stoma output (which includes digestive enzymes), or perspiration stay in contact with the skin.
- Dansac's **TRE technology** is designed to protect **and maintain** the skin's acid mantle, so that it can go on working to protect you.



► Colour Matching System

To help you select compatible 2 Piece pouches and wafers please use the colour codes to identify the matching ring sizes



36mm



55mm



43mm



70mm

Help keep your patients' skin at peace



With NEW NovaLife TRE skin barriers

Discover NovaLife TRE ostomy barriers from Dansac, specifically designed to help keep skin naturally healthy with three levels of protection.



What if a skin barrier could help keep skin naturally healthy?



Shelley, stoma since 2015

The best skin is healthy skin



NovaLife TRE ostomy skin barriers

Takes skin protection to another level. Three to be exact.

- Available in 1 piece closed, drainable and 2 piece wafers
- Stays in place, yet is easy to remove
- Designed to absorb stoma output and perspiration
- Maintaining normal skin pH helps protect it from damaging stoma output



To learn more or to obtain a sample, please contact your Dansac Territory Manager, Customer Engagement on **1800 880 840** or visit **www.dansac.com.au**

NovaLife TRE 1 Closed Flat



Max Issue 90 / month



Box Price \$90.27

NOVALIFE

	MIDI	MAXI	MAXI	MAXI
	EasiView™	EasiView™	EasiView™	Clear
	30 / box	30 / box	30 / box	30 / box
Starter hole	Code	Code	Code	Code
15 mm (CTF)	3001-15	3007-15*	3005-15	3006-15
25 mm (PC)	3001-25	3007-25*	3005-25	
30 mm (PC)	3001-30	3007-30*	3005-30	
35 mm (PC)	3001-35			
40 mm (PC)	3001-40			
45 mm (PC)	3001-45			

* MIDI barrier on MAXI pouch

Specifications

Max cut (h x w mm)	55 x 70	55 x 70	70 x 90	70 x 90
Pouch size (mm)	190 x 145	210 x 157	210 x 157	210 x 157
Volume capacity (ml)	375	555	530	530

NovaLife TRE 1 Closed Soft Convex



Max Issue 60 / month



Box Price \$47.07

NOVALIFE

		MIDI	MIDI	MAXI	MAXI
		EasiView™	EasiView™	EasiView™	EasiView™
		10 / box	10 / box	10 / box	10 / box
Starter hole	Plateau	Code	Code	Code	Code
15-24 mm (CTF)	A	3071-24			
15-34 mm (CTF)	B		3071-34		
15-44 mm (CTF)	C			3071-44	
15-54 mm (CTF)	D				3071-54
25 mm (PC)	A	3071-25			
30 mm (PC)	B		3071-30		
35 mm (PC)	B		3071-35		
40 mm (PC)	C			3071-40	

Specifications

Convexity depth (mm)	6	6	6	6
Skin barrier size (mm)	107 x 91	130 x 110	130 x 110	150 x 127
Pouch size (mm)	190 x 145	190 x 145	210 x 157	210 x 157
Volume capacity (ml)	450	430	590	570

NovaLife TRE 1 Closed Convex



Max Issue 60 / month



Box Price \$47.07

NOVALIFE

		MINI	MIDI	MIDI	MAXI
		Clear	EasiView™	EasiView™	EasiView™
		10 / box	10 / box	10 / box	10 / box
Starter hole	Plateau	Code	Code	Code	Code
15-24 mm (CTF)	A	3032-24			
15-37 mm (CTF)	B			3031-37	
15-59 mm (CTF)	D				3031-59
25 mm (PC)	A		3035-25		
30 mm (PC)	B			3031-30	
35 mm (PC)	B			3031-35	

Specifications

Convexity depth (mm)	6	6	6	6
Skin barrier size (mm)	107 x 91	107 x 91	116 x 98	150 x 127
Pouch size (mm)	175 x 130	190 x 145	190 x 145	210 x 157
Volume capacity (ml)	340	460	440	550

NovaLife TRE 1 Open Flat



Max Issue 30 / month

Box Price \$146.37 (MINI / MIDI)
Box Price \$48.79 (MAXI)

NOVALIFE

		MINI	MINI	MIDI	MIDI	MAXI	MAXI	MAXI
		Opaque	Clear	EasiView™	Clear	EasiView™	EasiView™	Clear
		30 / box	30 / box	30 / box	30 / box	10 / box	10 / box	10 / box
Starter hole		Code	Code	Code	Code	Code	Code	Code
15 mm (CTF)		3025-15	3026-15	3023-15	3024-15	3017-15*	3015-15	3016-15
25 mm (PC)				3023-25	3024-25	3017-25*	3015-25	
30 mm (PC)				3023-30	3024-30	3017-30*	3015-30	
35 mm (PC)				3023-35	3024-35			
40 mm (PC)				3023-40				

* MIDI barrier on MAXI pouch

Specifications

Max cut (h x w mm)	38 x 50	38 x 50	55 x 70	55 x 70	55 x 70	70 x 90	70 x 90
Skin barrier size (mm)	85 x 72	85 x 72	105 x 89	105 x 89	105 x 89	125 x 106	125 x 106
Pouch size (mm)	220 x 132	220 x 132	239 x 144	239 x 144	259 x 158	259 x 158	259 x 158
Volume capacity (ml)	245	245	345	345	535	445	445

NovaLife TRE 1 Open Soft Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

Starter hole	Plateau	MIDI	MIDI	MAXI	MAXI	MAXI	MAXI
		EasiView™	EasiView™	EasiView™	EasiView™	EasiView™	Clear
		10 / box	10 / box	10 / box	10 / box	10 / box	10 / box
		Code	Code	Code	Code	Code	Code
15-24 mm (CTF)	A	3081-24		3083-24			
15-34 mm (CTF)	B		3081-34	3083-34			3084-34
15-44 mm (CTF)	C				3081-44		
15-54 mm (CTF)	D					3081-54	
25 mm (PC)	A	3081-25		3083-25			
30 mm (PC)	B		3081-30	3083-30			
35 mm (PC)	B		3081-35	3083-35			
40 mm (PC)	C				3081-40		
45 mm (PC)	C				3081-45		

Specifications

Convexity depth (mm)	6	6	6	6	6	6
Skin barrier size (mm)	107 x 91	130 x 110	130 x 110	130 x 110	150 x 127	130 x 110
Pouch size (mm)	239 x 144	239 x 144	259 x 158	259 x 158	259 x 158	259 x 158
Volume capacity (ml)	350	360	540	520	500	540

NovaLife TRE 1 Open Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

Starter hole	Plateau	MIDI	MIDI	MIDI	MAXI
		EasiView™	EasiView™	EasiView™	EasiView™
		10 / box	10 / box	10 / box	10 / box
		Code	Code	Code	Code
15-37 mm (CTF)	B		3041-37		
15-46 mm (CTF)	C			3041-46	
15-59 mm (CTF)	D				3041-59
25 mm (PC)	A	3045-25			
30 mm (PC)	B		3041-30		
35 mm (PC)	B		3041-35		
40 mm (PC)	D			3041-40	

Specifications

Convexity depth (mm)	6	6	6	6
Skin barrier size (mm)	107 x 91	116 x 98	130 x 110	150 x 127
Pouch size (mm)	239 x 144	239 x 144	239 x 144	259 x 158
Volume capacity (ml)	420	380	350	500

NovaLife TRE 1 Open Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

		MINI	MIDI	MAXI
		Clear	Clear	Clear
		10 / box	10 / box	10 / box
Starter hole	Plateau	Code	Code	Code
15-24 mm (CTF)	A	3042-24		
15-37 mm (CTF)	B		3042-37	
15-59 mm (CTF)	D			3042-59
30 mm (PC)	B		3042-30	

Specifications

Convexity depth (mm)	6	6	6
Skin barrier size (mm)	107 x 91	116 x 98	150 x 127
Pouch size (mm)	220 x 132	239 x 144	259 x 158
Volume capacity (ml)	270	380	500

NovaLife TRE 2 Flat Wafer



Max Issue 30 / month



Box Price \$20.70

Refer to page 30-31 for corresponding NovaLife pouches

NOVALIFE

			5 / box	5 / box	5 / box	5 / box
Ring size	Starter hole	Max cut	Code	Code	Code	Code
36 mm	15 mm (CTF)	28 mm	3136-15			
43 mm	15 mm (CTF)	35 mm		3143-15		
	25 mm (PC)			3143-25		
	30 mm (PC)			3143-30		
55 mm	15 mm (CTF)	47 mm			3155-15	
	35 mm (PC)				3155-35	
	40 mm (PC)				3155-40	
70 mm	15 mm (CTF)	62 mm				3170-15

Specifications

Skin barrier size (mm)			85 x 72	105 x 89	105 x 89	125 x 106
------------------------	--	--	---------	----------	----------	-----------

NovaLife TRE 2 Convex Wafer



Max Issue 20 / month



Box Price \$29.20


Refer to page 30-31 for corresponding NovaLife pouches

NOVALIFE

Ring size	Plateau	Starter hole	Max cut	5 / box	5 / box	5 / box	5 / box
				Code	Code	Code	Code
36 mm	A	15 mm (CTF)	23 mm	3236-15			
43 mm	B	15 mm (CTF)	30 mm		3243-15		
		25 mm (PC)			3243-25		
		30 mm (PC)			3243-30		
55 mm	C	15 mm (CTF)	42 mm			3255-15	
		35 mm (PC)				3255-35	
		40 mm (PC)				3255-40	
70 mm	D	15 mm (CTF)	59 mm				3270-15

Specifications

Convexity depth (mm)			6	6	6	6
Skin barrier size (mm)			107 x 91	116 x 98	130 x 110	150 x 127



Dansac NovaLife
GX & GX+

Scott, ileostomy since 2006

NovaLife GX 1 Closed Flat



Max Issue 90 / month

Box Price \$90.27

NOVALIFE	MINI	MIDI	MIDI	MAXI	MAXI
	Opaque	EasiView™	Clear	EasiView™	Clear
	30 / box	30 / box	30 / box	30 / box	30 / box
Starter hole	Code	Code	Code	Code	Code
10 mm (CTF)	903-10	9010-10	902-10	9050-10	906-10
25 mm (PC)	903-25	9010-25	902-25	9050-25	
30 mm (PC)	903-30	9010-30	902-30	9050-30	
35 mm (PC)	903-35	9010-35	902-35	9050-35	
40 mm (PC)		9010-40	902-40		
45 mm (PC)		9010-45	902-45		
50 mm (PC)		9010-50	902-50		

Specifications

Max cut (h x w mm)	38 x 50	55 x 70	55 x 70	70 x 90	70 x 90
Skin barrier size (mm)	85 x 72	105 x 89	105 x 89	125 x 106	125 x 106
Pouch size (mm)	175 x 130	190 x 145	190 x 145	210 x 157	210 x 157
Volume capacity (ml)	360	470	470	620	620

Colour Matching System

To help you select compatible 2 Piece pouches and wafers please use the colour codes to identify the matching ring sizes

- 36mm
- 55mm
- 43mm
- 70mm

NovaLife GX 1 Closed Soft Convex



Max Issue 60 / month



Box Price \$47.07

NOVALIFE

		MIDI	MIDI	MAXI	MAXI
		EasiView™	EasiView™	EasiView™	EasiView™
		10 / box	10 / box	10 / box	10 / box
Starter hole	Plateau	Code	Code	Code	Code
15-24 mm (CTF)	A	971-24			
15-34 mm (CTF)	B		971-34		
15-44 mm (CTF)	C			971-44	
15-54 mm (CTF)	D				971-54
25 mm (PC)	A	971-25			
30 mm (PC)	B		971-30		
35 mm (PC)	B		971-35		
40 mm (PC)	C			971-40	

Specifications

Convexity depth (mm)	6	6	6	6
Skin barrier size (mm)	107 x 91	130 x 110	130 x 110	150 x 127
Pouch size (mm)	190 x 145	190 x 145	210 x 157	210 x 157
Volume capacity (ml)	450	430	590	570

NovaLife GX 1 Closed Convex



Max Issue 60 / month



Box Price \$47.07

NOVALIFE

		MIDI	MIDI	MIDI	MIDI	MAXI
		EasiView™	EasiView™	Clear	Clear	EasiView™
		10 / box	10 / box	10 / box	10 / box	10 / box
Starter hole	Plateau	Code	Code	Code	Code	Code
25 mm (PC)	A	9310-25				
30 mm (CTF)	B		9310-30	932-37		
35 mm (PC)	B		9310-35	932-35		
45 mm (PC)	C				932-45	
50 mm (CTF)	D					9310-50

Specifications

Convexity depth (mm)	6	6	6	6	6
Skin barrier size (mm)	107 x 91	116 x 98	116 x 98	130 x 110	150 x 127
Pouch size (mm)	190 x 145	190 x 145	190 x 145	190 x 145	210 x 157
Volume capacity (ml)	460	440	440	430	550

NovaLife GX 1 Open Flat



Max Issue 30 / month



Box Price \$146.37 (MINI / MIDI)
Box Price \$48.79 (MAXI)

NOVALIFE

	MINI	MINI	MIDI	MIDI	MAXI	MAXI
	Opaque	Clear	Opaque	Clear	Opaque	Clear
	30 / box	30 / box	30 / box	30 / box	10 / box	10 / box
Starter hole	Code	Code	Code	Code	Code	Code
10 mm (CTF)	925-10	926-10	923-10	924-10	915-10	916-10
25 mm (PC)			923-25	924-25	915-25	
30 mm (PC)			923-30	924-30	915-30	916-30
35 mm (PC)			923-35	924-35	915-35	916-35
40 mm (PC)			923-40	924-40		
45 mm (PC)			923-45	924-45		

Specifications

Max cut (h x w mm)	38 x 50	38 x 50	55 x 70	55 x 70	70 x 90	70 x 90
Skin barrier size (mm)	85 x 72	85 x 72	105 x 89	105 x 89	125 x 106	125 x 106
Pouch size (mm)	220 x 132	220 x 132	239 x 144	239 x 144	259 x 158	259 x 158
Volume capacity (ml)	300	300	420	420	550	550

NovaLife GX+ 1 Open Flat



Max Issue 30 / month



Box Price \$146.37 (MIDI)
Box Price \$48.79 (MAXI)

NOVALIFE

	MIDI	MAXI	MAXI
	EasiView™	EasiView™	Clear
	30 / box	10 / box	10 / box
Starter hole	Code	Code	Code
10 mm (CTF)	9630-10	9610-10	962-10
25 mm (PC)	9630-25	9610-25	
30 mm (PC)	9630-30	9610-30	
35 mm (PC)		9610-35	

Specifications

Max cut (h x w mm)	55 x 70	70 x 90	70 x 90
Skin barrier size (mm)	105 x 89	125 x 106	125 x 106
Pouch size (mm)	239 x 144	259 x 158	259 x 158
Volume capacity (ml)	420	550	550

NovaLife GX 1 Open Soft Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

MIDI

MIDI

EasiView™

EasiView™

10 / box

10 / box

Starter hole

Plateau

Code

Code

15-24 mm (CTF)	A	981-24	
15-34 mm (CTF)	B		981-34
15-44 mm (CTF)	C		
15-54 mm (CTF)	D		
25 mm (PC)	A	981-25	
30 mm (PC)	B		981-30
35 mm (PC)	B		981-35
40 mm (PC)	C		
45 mm (PC)	C		

Specifications

Convexity depth (mm)	6	6
Skin barrier size (mm)	107 x 91	130 x 110
Pouch size (mm)	239 x 144	239 x 144
Volume capacity (ml)	350	360

NovaLife GX 1 Open Soft Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

MAXI

MAXI

MAXI

MAXI

MAXI

Clear

Clear

EasiView™

EasiView™

EasiView™

10 / box

10 / box

10 / box

10 / box

10 / box

Starter hole

Plateau

Code

Code

Code

Code

Code

15-24 mm (CTF)	A	984-24		983-24		
15-34 mm (CTF)	B	984-34		983-34		
15-44 mm (CTF)	C		982-44		981-44	
15-54 mm (CTF)	D					981-54
25 mm (PC)	A	984-25		983-25		
30 mm (PC)	B	984-30		983-30		
35 mm (PC)	B	984-35		983-35		
40 mm (PC)	C		982-40		981-40	
45 mm (PC)	C		982-45		981-45	

Specifications

Convexity depth (mm)	6	6	6	6	6
Skin barrier size (mm)	130 x 110	130 x 110	130 x 110	130 x 110	150 x 127
Pouch size (mm)	259 x 158	259 x 158	259 x 158	259 x 158	259 x 158
Volume capacity (ml)	540	520	540	520	500

NovaLife GX 1 Open Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

Starter hole	Plateau	MINI	MIDI	MIDI	MIDI	MIDI	MAXI	MAXI
		Opaque	Opaque	Clear	Opaque	Clear	Opaque	Clear
		10 / box	10 / box	10 / box	10 / box	10 / box	10 / box	10 / box
		Code	Code	Code	Code	Code	Code	Code
10-37 mm (CTF)	B		941-37	942-37				
10-46 mm (CTF)	C				941-46			
10-59 mm (CTF)	D						941-59	942-59
20 mm (PC)	A	941-20						
25 mm (PC)	A	941-25						
30 mm (PC)	B		941-30	942-30				
35 mm (PC)	B		941-35	942-35				
40 mm (PC)	C				941-40	942-40		
45 mm (PC)	C					942-45		
50 mm (PC)	D						941-50	942-50

Specifications

Convexity depth (mm)	6	6	6	6	6	6	6
Skin barrier size (mm)	107 x 91	116 x 98	116 x 98	130 x 110	130 x 110	150 x 127	150 x 127
Pouch size (mm)	220 x 132	239 x 144	239 x 144	239 x 144	239 x 144	259 x 158	259 x 158
Volume capacity (ml)	270	380	380	350	350	500	500

NovaLife GX+ 1 Open Convex



Max Issue 30 / month



Box Price \$65.78

NOVALIFE

Starter hole	Plateau	MIDI	MIDI	MAXI	MAXI
		EasiView™	EasiView™	EasiView™	Clear
		10 / box	10 / box	10 / box	10 / box
		Code	Code	Code	Code
10-59 mm (CTF)	D			9510-59	952-59
20 mm (PC)	A	9510-20			
25 mm (PC)	A	9510-25			
30 mm (PC)	B		9510-30		
35 mm (PC)	B		9510-35		
50 mm (PC)	D			9510-50	

Specifications

Convexity depth (mm)	6	6	6	6
Skin barrier size (mm)	107 x 91	116 x 98	150 x 127	150 x 127
Pouch size (mm)	239 x 144	239 x 144	259 x 158	259 x 158
Volume capacity (ml)	420	380	500	500

NovaLife 2 Closed Pouch



Max Issue 90 / month



Box Price \$53.07

Refer to page 17-18 for NovaLife TRE 2 piece wafers

NOVALIFE

	MINI	MIDI	MIDI
	Opaque	Opaque	Clear
	30 / box	30 / box	30 / box
Ring size	Code	Code	Code
36 mm	1303-36	1301-36	
43 mm	1303-43	1301-43	1302-43
55 mm	1303-55	1301-55	1302-55
70 mm		1301-70	1302-70

Specifications

Pouch size (mm)	175 x 130	190 x 145	190 x 145
Volume capacity (ml)	260	390	390

NovaLife 2 Open Pouch



Max Issue 60 / month



Box Price \$111.27

Refer to page 17-18 for NovaLife TRE 2 piece wafers

NOVALIFE

	MINI	MINI	MIDI	MIDI	MAXI	MAXI
	Opaque	Clear	Opaque	Clear	Opaque	Clear
	30 / box	30 / box	30 / box	30 / box	30 / box	30 / box
Ring size	Code	Code	Code	Code	Code	Code
36 mm	1319-36	1320-36	1315-36	1316-36		
43 mm	1319-43		1315-43	1316-43	1321-43	1322-43
55 mm	1319-55		1315-55	1316-55	1321-55	1322-55
70 mm			1315-70	1316-70	1321-70	1322-70

Specifications

Pouch size (mm)	220 x 132	220 x 132	239 x 144	239 x 144	259 x 158	259 x 158
Volume capacity (ml)	280	300	440	480	650	690

NovaLife GX 2 Flat Wafer



Max Issue 30 / month



Box Price \$20.70

NOVALIFE

Ring size	Starter hole	Max cut	5 / box	5 / box	5 / box	5 / box
			Code	Code	Code	Code
36 mm	10 mm (CTF) 25 mm (PC)	28 mm	1736-10 1736-25			
43 mm	10 mm (CTF) 25 mm (PC) 30 mm (PC)	35 mm		1743-10 1743-25 1743-30		
55 mm	10 mm (CTF) 35 mm (PC) 40 mm (PC)	47 mm			1755-10 1755-35 1755-40	
70 mm	10 mm (CTF)	62 mm				1770-10

Specifications

Skin barrier size (mm)			85 x 72	105 x 89	105 x 89	125 x 106
------------------------	--	--	---------	----------	----------	-----------

NovaLife GX+ 2 Flat Wafer



Max Issue 20 / month



Box Price \$23.71

NOVALIFE

Ring size	Starter hole	Max cut	5 / box	5 / box	5 / box
			Code	Code	Code
43 mm	10 mm (CTF) 25 mm (PC)	35 mm	1843-10 1843-25		
55 mm	10 mm (CTF)	47 mm		1855-10	
70 mm	10 mm (CTF) 45 mm (PC)	62 mm			1870-10 1870-45

Specifications

Skin barrier size (mm)			105 x 89	105 x 89	125 x 106
------------------------	--	--	----------	----------	-----------

NovaLife GX 2 Convex Wafer



Max Issue 20 / month



Box Price \$29.20

NOVALIFE

Ring size	Plateau	Starter hole	Max cut	5 / box			
				Code	Code	Code	Code
36 mm	A	10 mm (CTF) 20 mm (PC) 23 mm (PC)	23 mm	1936-10 1936-20 1936-23			
43 mm	B	10 mm (CTF) 25 mm (PC) 30 mm (PC)	30 mm		1943-10 1943-25 1943-30		
55 mm	C	10 mm (CTF) 35 mm (PC) 40 mm (PC)	42 mm			1955-10 1955-35 1955-40	
70 mm	D	10 mm (CTF) 40 mm (PC) 45 mm (PC)	59 mm				1970-10 1970-40 1970-45

Specifications

Convexity depth (mm)			6	6	6	6
Skin barrier size (mm)			107 x 91	116 x 98	130 x 110	150 x 127

NovaLife GX+ 2 Convex Wafer



Max Issue 20 / month



Box Price \$29.20

NOVALIFE

Ring size	Plateau	Starter hole	Max cut	5 / box			
				Code	Code	Code	Code
36 mm	A	20 mm (PC)	23 mm	2036-20			
43 mm	B	10 mm (CTF) 25 mm (PC) 30 mm (PC)	30 mm		2043-10 2043-25 2043-30		
55 mm	C	10 mm (CTF) 35 mm (PC)	42 mm			2055-10 2055-35	
70 mm	D	10 mm (CTF) 45 mm (PC)	59 mm				2070-10 2070-45

Specifications

Convexity depth (mm)			6	6	6	6
Skin barrier size (mm)			107 x 91	116 x 98	130 x 110	150 x 127



Dansac Nova

Nova 1 Closed Flat



Max Issue 60 / month
Max Issue 90 / month

Box Price \$41.19 (MINI CAP)
Box Price \$82.05

NOVA¹

	MINI CAP	MINI	MIDI	MIDI	MAXI	MAXI
	Opaque	Opaque	Opaque	Clear	Opaque	Clear
	30 / box	30 / box	30 / box	30 / box	30 / box	30 / box
Starter hole	Code	Code	Code	Code	Code	Code
20-50 mm (CTF)	829-20	803-20				
20-60 mm (CTF)			801-20	802-20	805-20	806-20
25 mm (PC)	829-25	803-25	801-25	802-25		
30 mm (PC)	829-30	803-30	801-30	802-30		
35 mm (PC)		803-35	801-35	802-35		
40 mm (PC)	829-40	803-40	801-40	802-40		
45 mm (PC)			801-45	802-45		
50 mm (PC)	829-50		801-50	802-50		

Specifications

Skin barrier size (mm)	92	85 x 85	100 x 100	100 x 100	100 x 100	100 x 100
Pouch size (mm)	105	164 x 128	194 x 150	194 x 150	217 x 150	217 x 150
Volume capacity (ml)	50	230	420	420	560	560

Colour Matching System

To help you select compatible 2 Piece pouches and wafers please use the colour codes to identify the matching ring sizes

- 36mm
- 55mm
- 43mm
- 70mm

Nova 1 Closed Soft Mouldable



- 5mm soft convexity depth designed to give gentle pressure
- GX hydrocolloid barrier is designed for skin-friendliness and flexibility
- Nova filter absorbs odours and is designed to reduce the risk of ballooning
- Soft and water-repellent backing



Max Issue 60 / month



Box Price \$44.33

NOVA¹

MIDI

Opaque

10 / box

Starter hole	Max stretch	Plateau	Code
25 mm	35 mm	A	871-25
30 mm	40 mm	B	871-30
35 mm	45 mm	B	871-35

Specifications

Convexity (mm)	5
Skin barrier size (mm)	106 x 106
Pouch size (mm)	194 x 150
Volume capacity (ml)	470

Nova 1 Closed Convex



Max Issue 60 / month



Box Price \$44.33

NOVA¹

MIDI

Opaque

10 / box

MAXI

Opaque

10 / box

MAXI

Clear

10 / box

Starter hole	Plateau	Code	Code	Code
15-37 mm (CTF)	B		831-37	832-37
15-46 mm (CTF)	C			832-46
25 mm (PC)	A	831-25		
30 mm (PC)	B		831-30	
35 mm (PC)	B		831-35	
40 mm (PC)	C		831-40	
45 mm (PC)	C		831-45	

Specifications

Convexity (mm)	6	6	6
Skin barrier size (mm)	100 x 100	115 x 115	115 x 115
Pouch size (mm)	194 x 150	217 x 150	217 x 150
Volume capacity (ml)	480	550	550

Nova 2 Closed Pouch



Max Issue 90 / month



Box Price \$46.92

NOVA²

	MINI	MIDI	MIDI
	Opaque	Opaque	Clear
	30 / box	30 / box	30 / box
Ring size	Code	Code	Code
36 mm		1201-36	
43 mm	1203-43	1201-43	1202-43
55 mm	1203-55	1201-55	1202-55
70 mm		1201-70	

Specifications

Pouch size (mm)	164 x 128	194 x 150	194 x 150
Volume capacity (ml)	200	450	450

Nova 1 FoldUp Flat



Max Issue 30 / month

Box Price \$140.22 (MINI / MIDI)
Box Price \$46.74 (MAXI)

NOVA¹ FoldUp

	MINI	MINI	MIDI	MIDI	MAXI	MAXI
	Opaque	Clear	Opaque	Clear	Opaque	Clear
	30 / box	30 / box	30 / box	30 / box	10 / box	10 / box
Starter hole	Code	Code	Code	Code	Code	Code
15-50 mm (CTF)	825-15	826-15*				
15-60 mm (CTF)			823-15	824-15*		
15-90 mm (CTF)					815-15	816-15*
20 mm (PC)			823-20			
25 mm (PC)			823-25			
30 mm (PC)			823-30			
35 mm (PC)			823-35			
40 mm (PC)			823-40			
45 mm (PC)			823-45			

* Post op pouches

Specifications

Max cut (h x w mm)	50 x 50	50 x 50	60 x 60	60 x 60	90 x 90	90 x 90
Skin barrier size (mm)	85 x 85	85 x 85	100 x 100	100 x 100	125 x 125	125 x 125
Pouch size (mm)	225 x 128	225 x 128	265 x 148	265 x 148	306 x 180	306 x 180
Volume capacity (ml)	290	290	510	510	800	800

Nova 1 FoldUp X3



Max Issue 30 / month



Box Price \$63.73

NOVA¹
FoldUp

	MIDI	MIDI
	Opaque	Clear
	10 / box	10 / box
Starter hole	Code	Code
15-46 mm (CTF)	861-15	862-15
25 mm (PC)	861-25	
30 mm (PC)	861-30	
35 mm (PC)	861-35	

Specifications

Skin barrier size (mm)	106 x 106	106 x 106
Pouch size (mm)	265 x 148	265 x 148
Volume capacity (ml)	540	540

Nova 1 FoldUp Soft Mouldable



Max Issue 30 / month



Box Price \$63.73

NOVA¹

			MIDI	MIDI
			Opaque	Clear
			10 / box	10 / box
Starter hole	Max stretch	Plateau	Code	Code
25 mm	35 mm	A	881-25	882-25
30 mm	40 mm	B	881-30	882-30
35 mm	45 mm	B	881-35	882-35

Specifications

Convexity (mm)	5	5
Skin barrier size (mm)	106 x 106	106 x 106
Pouch size (mm)	265 x 148	265 x 148
Volume capacity (ml)	550	550

Nova 1 FoldUp Convex



Max Issue 30 / month



Box Price \$63.73

NOVA¹
FoldUp

		MIDI		MIDI	
		Opaque	Clear	Opaque	Clear
		10 / box	10 / box	10 / box	10 / box
		Code	Code	Code	Code
Starter hole	Plateau				
15-24 mm (CTF)	A	841-24	842-24		
15-37 mm (CTF)	B	841-37	842-37		
15-46 mm (CTF)	C			841-46	842-46
20 mm (PC)	A	841-20			
25 mm (PC)	A	841-25	842-25		
30 mm (PC)	B	841-30	842-30		
35 mm (PC)	B	841-35	842-35		

Specifications

Convexity (mm)	6	6	6	6
Skin barrier size (mm)	100 x 100	100 x 100	115 x 115	115 x 115
Pouch size (mm)	265 x 148	265 x 148	265 x 148	265 x 148
Volume capacity (ml)	550	550	550	550

Nova 1 Infant / Paediatric



Max Issue 60 / month



Box Price \$100.47

NOVA¹
Infant

	INFANT	
	Clear	Clear
	30 / box	30 / box
	Code	Code
Starter hole		
0-40 mm (CTF)*	818-00	
10-40 mm (CTF)		818-10

* No starter hole

Specifications

Skin barrier size (mm)	70 x 70	70 x 70
Pouch size (mm)	208 x 111	208 x 111
Volume capacity (ml)	250	250
Soft Wire Ties	30	30

Nova 2 FoldUp Pouch



Max Issue 60 / month



Box Price \$35.04

NOVA²
FoldUp

	MINI	MIDI	MIDI	MAXI	MAXI
	Opaque	Opaque	Clear	Opaque	Clear
	10 / box	10 / box	10 / box	10 / box	10 / box
Ring size	Code	Code	Code	Code	Code
36 mm	1219-36	1215-36	1216-36		
43 mm	1219-43	1215-43	1216-43		
55 mm	1219-55	1215-55	1216-55	1221-55	1222-55
70 mm		1215-70	1216-70	1221-70	1222-70

Specifications

Pouch size (mm)	225 x 128	265 x 148	265 x 148	294 x 155	294 x 155
Volume capacity (ml)	330	590	590	820	820

Nova 2 High Output Pouch



Max Issue 20 / month



Box Price \$163.80

NOVA²
HighOutput

	MIDI	MIDI
	Opaque	Clear
	10 / box	10 / box
Ring size	Code	Code
43 mm	1207-43	1208-43
55 mm	1207-55	1208-55
70 mm	1207-70	1208-70

Specifications

Pouch size (mm)	294 x 155	294 x 155
Volume capacity (ml)	840	840

Nova 1 Urostomy X3



Max Issue 30 / month

Box Price \$63.73 (Pouch)
Box Price \$39.12 (Adaptor)

NOVA¹ X3

Urostomy

	POUCH	POUCH	ADAPTOR
	Opaque	Clear	Clear
	10 / box*	10 / box*	10 / box*
Starter hole	Code	Code	Code
			400-50
12-46 mm (CTF)	891-12	892-12	
20 mm (PC)	891-20		
25 mm (PC)	891-25		
30 mm (PC)	891-30		

Specifications

Skin barrier size (mm)	106 x 106	106 x 106	
Pouch size (mm)	245 x 150	245 x 150	
Volume capacity (ml)	330	330	

* + 1 adaptor and 10 caps

Nova 1 Urostomy Convex



Max Issue 30 / month

Box Price \$63.73 (Pouch)
Box Price \$39.12 (Adaptor)

NOVA¹

Urostomy

		POUCH	POUCH	ADAPTOR
		Clear	Clear	Clear
		10 / box*	10 / box*	10 / box*
Starter hole	Plateau	Code	Code	Code
				400-50
15-24 mm (CTF)	A	894-24		
15-37 mm (CTF)	B	894-37		
15-46 mm (CTF)	C		894-46	
20 mm	A	894-20		
25 mm	A	894-25		
30 mm	B	894-30		

Specifications

Convexity (mm)	6	6	
Skin barrier size (mm)	100 x 100	115 x 115	
Pouch size (mm)	245 x 150	245 x 150	
Volume capacity (ml)	330	330	

* + 1 adaptor and 10 caps

Nova 2 Urostomy Pouch



Max Issue 30 / month



Box Price \$35.04 (Pouch)
Box Price \$39.12 (Adaptor)

NOVA²
Urostomy

	POUCH	POUCH	ADAPTOR
	Opaque	Clear	Clear
	10 / box*	10 / box*	10 / box*
Starter hole	Code	Code	Code
			400-50
 36 mm	1217-36	1218-36	
 43 mm	1217-43	1218-43	
 55 mm	1217-55	1218-55	

Specifications

Pouch size (mm)	245 x 150	245 x 150	
Volume capacity (ml)	330	330	

* + 1 adaptor and 10 caps

Nova 2 Flat Wafer



Max Issue 30 / month



Box Price \$20.70

NOVA²
Wafer

			5 / box
Ring size	Starter hole	Max cut	Code
 36 mm 85 x 85 mm	20 mm (PC) 25 mm (PC)	28 mm	1136-20 1136-25
 43 mm 100 x 100 mm	15 mm (CTF) 25 mm (PC) 30 mm (PC)	35 mm	1143-15 1143-25 1143-30
 55 mm 100 x 100 mm	10 mm (CTF) 35 mm (PC) 40 mm (PC)	47 mm	1155-15 1155-35 1155-40
 70 mm 115 x 115 mm	15 mm (CTF) 45 mm (PC) 50 mm (PC)	62 mm	1170-15 1170-45 1170-50

Nova 2 X3 Wafer



 Max Issue 20 / month

 Box Price \$29.20

NOVA² X3

Wafer

			5 / box
Ring size	Starter hole	Max cut	Code
 43 mm 106 x 106 mm	15 mm (CTF)	30 mm	1343-15
	25 mm (PC)	25 mm	1343-25
 55 mm 121 x 121 mm	15 mm (CTF)	42 mm	1355-15
	32 mm (PC)	40 mm	1355-32
	35 mm (PC)	40 mm	1355-35

Nova 2 Soft Mouldable Wafer



 Max Issue 20 / month

 Box Price \$29.20

NOVA²

Soft convex Wafer

				5 / box
Ring size	Plateau	Starter hole	Max stretch	Code
 43 mm 106 x 106 mm	C	18 mm (PC)	28 mm	1443-18
		21 mm (PC)	31 mm	1443-21
		25 mm (PC)	35 mm	1443-25
		28 mm (PC)	38 mm	1443-28
 55 mm 121 x 121 mm	D	32 mm (PC)	42 mm	1455-32
		35 mm (PC)	45 mm	1455-35
		38 mm (PC)	48 mm	1455-38

Specifications

Convexity (mm)			5
----------------	--	--	---

Nova 2 Convex Wafer



Max Issue 20 / month




Box Price \$29.20

NOVA²
Convex W a f e r

				5 / box
Ring size	Plateau	Starter hole	Max cut	Code
36 mm 100 x 100 mm	B	15 mm (CTF)	23 mm	1536-15
		20 mm (PC)	20 mm	1536-20
		23 mm (PC)	23 mm	1536-23
43 mm 100 x 100 mm	C	15 mm (CTF)	30 mm	1543-15
		25 mm (PC)	25 mm	1543-25
		30 mm (PC)	30 mm	1543-30
55 mm 115 x 115 mm	D	15 mm (CTF)	42 mm	1555-15
		35 mm (PC)	35 mm	1555-35
		40 mm (PC)	40 mm	1555-40

Specifications

Convexity (mm)			6
----------------	--	--	---



Dansac Accessories

Mariam, colostomy since 2006

Dansac TRE Seals

 Max Issue 30 / month

 Box Price \$45.31 (072-48)
Box Price \$135.93 (071-30 & 071-40)

Starter hole	Diameter size	Code	Units
18 mm	48 mm	072-48	10 / box
30 mm	57.5 mm	071-30	30 / box
40 mm	70 mm	071-40	30 / box



Dansac Seals

 Max Issue 30 / month

 Box Price \$135.93

Starter hole	Code	Units
20 mm	070-20	30 / box
30 mm	070-30	30 / box
40 mm	070-40	30 / box
50 mm	070-50	30 / box



Dansac GX-tra Seals

 Max Issue 30 / month

 Box Price \$135.93

Starter hole	Code	Units
20 mm	725-20/30	30 / box
30 mm	725-30/30	30 / box
40 mm	725-40/30	30 / box
50 mm	725-50/30	30 / box



Dansac X-tra Strips

Skin-friendly hydrocolloid barrier extender

 Max Issue 60 / month  Box Price \$43.74

Size	Code	Units
56.5 x 128 mm	075-30	30 / box



Dansac EasiSpray™ Adhesive Remover

Quick and easy removal of appliance

 Max Issue 2 / month  Box Price \$10.23

Volume	Code
50 ml	083-01



Dansac Soft Paste

Levelling skin folds and scars

 Max Issue 3 / month  Box Price \$11.28

Volume	Code
50 g	77550-0



Dansac Skin Lotion Tissue

Efficient cleansing of peristomal skin

 Max Issue 2 / month
Max Issue 5 / month  Box Price \$23.90 (Wipes)
Box Price \$4.77 (Travel Pack)

	Code	Units
Individual Wipes	71000-0000	50 / box
Travel Pack	71000-0010	10 / pack




Dansac Skin Lotion

Efficient cleansing of peristomal skin

 Max Issue 1 / month  Box Price \$8.35

Volume	Code
200 ml	70000-0000



Dansac Skin Creme

Extra moisture for your skin

 Max Issue 2 / month  Box Price \$7.92

Volume	Code
100 ml	085-00



Dansac Nodor S

Odour prevention for extra discretion

 Max Issue 2 / month  Box Price \$4.05

Volume	Code
50 ml	080-00



Belts

 Max Issue 4 / annum  Box Price \$5.98

	BELT 100 CM	BELT 150 CM
	OPAQUE	OPAQUE
Box	Code	Code
Box of 1 pcs	090-02	090-03



Clamps

 Max Issue 10 / annum  Box Price \$57.50

Box	Code
Box of 25 pcs	092-00



Dansac Irrigation Set

 900-20 Irrigation set – with ring holder and silicone ring (for 2-piece systems)

 900-35 Irrigation set – with hydrocolloid skin barrier

Content	Code	Price	Max Issue	Irrigation Set No.	
Water container	95200-0000	\$25.26	2 / annum		
Cone with tubing	95205-0000	\$12.97	1 / annum		
Dansac Skin Lotion 50 ml	70005-0001				
Dansac Skin Lotion Tissues	71000-0000	\$23.90	2 / month		
Dansac NovaLife 1 Mini 10 mm (CTF)	903-10				
Dansac Nova 1 Mini Cap 30 mm (PC)	829-20				
Irrigation sleeve with ring holder and silicone ring	950-20	\$26.18	50 / month		
Irrigation sleeve with hydrocolloid skin barrier	950-35	\$26.18	50 / month		
Belt 100 cm (medium)	090-02	\$5.98	4 / annum		

NEW

SEALING IN SKIN HEALTH



Dansac TRE Seal is more than just a seal – with three levels of protection, the Dansac TRE Seal has been designed to help keep healthy skin naturally healthy.



Adhesion

Designed to provide a secure, flexible seal to protect the skin from stoma fluid and to be easy to remove.



Absorption

Helps absorb excess moisture to maintain skin's natural balance without the seal losing internal or external strength.



pH balance

Designed to help manage the skin-damaging effects of digestive enzyme activity.

To learn more or to obtain a sample, please contact your Dansac Territory Manager, Customer Engagement on **1800 880 840** or visit **www.dansac.com.au**

Index

Index

NovaLife TRE

Code	SAS Code	Description	Page
3001-15	80194R	NovaLife TRE 1 Closed Flat	10
3001-25	80194R	NovaLife TRE 1 Closed Flat	10
3001-30	80194R	NovaLife TRE 1 Closed Flat	10
3001-35	80194R	NovaLife TRE 1 Closed Flat	10
3001-40	80194R	NovaLife TRE 1 Closed Flat	10
3001-45	80194R	NovaLife TRE 1 Closed Flat	10
3005-15	80194R	NovaLife TRE 1 Closed Flat	10
3005-25	80194R	NovaLife TRE 1 Closed Flat	10
3005-30	80194R	NovaLife TRE 1 Closed Flat	10
3006-15	80194R	NovaLife TRE 1 Closed Flat	10
3007-15	80194R	NovaLife TRE 1 Closed Flat	10
3007-25	80194R	NovaLife TRE 1 Closed Flat	10
3007-30	80194R	NovaLife TRE 1 Closed Flat	10
3015-15	80192P	NovaLife TRE 1 Open Flat	13
3015-25	80192P	NovaLife TRE 1 Open Flat	13
3015-30	80192P	NovaLife TRE 1 Open Flat	13
3016-15	80192P	NovaLife TRE 1 Open Flat	13
3017-15	80192P	NovaLife TRE 1 Open Flat	13
3017-25	80192P	NovaLife TRE 1 Open Flat	13
3017-30	80192P	NovaLife TRE 1 Open Flat	13
3023-15	80195T	NovaLife TRE 1 Open Flat	13
3023-25	80195T	NovaLife TRE 1 Open Flat	13
3023-30	80195T	NovaLife TRE 1 Open Flat	13
3023-35	80195T	NovaLife TRE 1 Open Flat	13
3023-40	80195T	NovaLife TRE 1 Open Flat	13
3024-15	80195T	NovaLife TRE 1 Open Flat	13
3024-25	80195T	NovaLife TRE 1 Open Flat	13
3024-30	80195T	NovaLife TRE 1 Open Flat	13
3024-35	80195T	NovaLife TRE 1 Open Flat	13
3025-15	80195T	NovaLife TRE 1 Open Flat	13
3026-15	80195T	NovaLife TRE 1 Open Flat	13
3031-30	80200C	NovaLife TRE 1 Closed Convex	12
3031-35	80200C	NovaLife TRE 1 Closed Convex	12
3031-37	80200C	NovaLife TRE 1 Closed Convex	12
3031-59	80200C	NovaLife TRE 1 Closed Convex	12
3032-24	80200C	NovaLife TRE 1 Closed Convex	12
3035-25	80200C	NovaLife TRE 1 Closed Convex	12
3041-30	80199B	NovaLife TRE 1 Open Convex	15
3041-35	80199B	NovaLife TRE 1 Open Convex	15
3041-37	80199B	NovaLife TRE 1 Open Convex	15
3041-40	80199B	NovaLife TRE 1 Open Convex	15
3041-46	80199B	NovaLife TRE 1 Open Convex	15
3041-59	80199B	NovaLife TRE 1 Open Convex	15
3042-24	80199B	NovaLife TRE 1 Open Convex	16

Index

NovaLife TRE

Code	SAS Code	Description	Page
3042-30	80199B	NovaLife TRE 1 Open Convex	16
3042-37	80199B	NovaLife TRE 1 Open Convex	16
3042-59	80199B	NovaLife TRE 1 Open Convex	16
3045-25	80199B	NovaLife TRE 1 Open Convex	15
3071-24	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-25	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-30	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-34	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-35	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-40	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-44	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3071-54	80193Q	NovaLife TRE 1 Closed Soft Convex.....	11
3081-24	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-25	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-30	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-34	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-35	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-40	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-44	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-45	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3081-54	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3083-24	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3083-25	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3083-30	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3083-34	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3083-35	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3084-34	80191N	NovaLife TRE 1 Open Soft Convex.....	14
3136-15	80201D	NovaLife TRE 2 Flat Wafer.....	17
3143-15	80201D	NovaLife TRE 2 Flat Wafer.....	17
3143-25	80201D	NovaLife TRE 2 Flat Wafer.....	17
3143-30	80201D	NovaLife TRE 2 Flat Wafer.....	17
3155-15	80201D	NovaLife TRE 2 Flat Wafer.....	17
3155-35	80201D	NovaLife TRE 2 Flat Wafer.....	17
3155-40	80201D	NovaLife TRE 2 Flat Wafer.....	17
3170-15	80201D	NovaLife TRE 2 Flat Wafer.....	17
3236-15	80202E	NovaLife TRE 2 Convex Wafer.....	18
3243-15	80202E	NovaLife TRE 2 Convex Wafer.....	18
3243-25	80202E	NovaLife TRE 2 Convex Wafer.....	18
3243-30	80202E	NovaLife TRE 2 Convex Wafer.....	18
3255-15	80202E	NovaLife TRE 2 Convex Wafer.....	18
3255-35	80202E	NovaLife TRE 2 Convex Wafer.....	18
3255-40	80202E	NovaLife TRE 2 Convex Wafer.....	18
3270-15	80202E	NovaLife TRE 2 Convex Wafer.....	18

Index

NovaLife GX & GX+

Code	SAS Code	Description	Page
902-10	5672M	NovaLife GX 1 Closed Flat	21
902-25	5672M	NovaLife GX 1 Closed Flat	21
902-30	5672M	NovaLife GX 1 Closed Flat	21
902-35	5672M	NovaLife GX 1 Closed Flat	21
902-40	5672M	NovaLife GX 1 Closed Flat	21
902-45	5672M	NovaLife GX 1 Closed Flat	21
902-50	5672M	NovaLife GX 1 Closed Flat	21
903-10	5672M	NovaLife GX 1 Closed Flat	21
903-25	5672M	NovaLife GX 1 Closed Flat	21
903-30	5672M	NovaLife GX 1 Closed Flat	21
903-35	5672M	NovaLife GX 1 Closed Flat	21
906-10	5672M	NovaLife GX 1 Closed Flat	21
915-10	3954D	NovaLife GX 1 Open Flat.....	24
915-25	3954D	NovaLife GX 1 Open Flat.....	24
915-30	3954D	NovaLife GX 1 Open Flat.....	24
915-35	3954D	NovaLife GX 1 Open Flat.....	24
916-10	3954D	NovaLife GX 1 Open Flat.....	24
916-30	3954D	NovaLife GX 1 Open Flat.....	24
916-35	3954D	NovaLife GX 1 Open Flat.....	24
923-10	3931X	NovaLife GX 1 Open Flat.....	24
923-25	3931X	NovaLife GX 1 Open Flat.....	24
923-30	3931X	NovaLife GX 1 Open Flat.....	24
923-35	3931X	NovaLife GX 1 Open Flat.....	24
923-40	3931X	NovaLife GX 1 Open Flat.....	24
923-45	3931X	NovaLife GX 1 Open Flat.....	24
924-10	3931X	NovaLife GX 1 Open Flat.....	24
924-25	3931X	NovaLife GX 1 Open Flat.....	24
924-30	3931X	NovaLife GX 1 Open Flat.....	24
924-35	3931X	NovaLife GX 1 Open Flat.....	24
924-40	3931X	NovaLife GX 1 Open Flat.....	24
924-45	3931X	NovaLife GX 1 Open Flat.....	24
925-10	3931X	NovaLife GX 1 Open Flat.....	24
926-10	3931X	NovaLife GX 1 Open Flat.....	24
932-35	3930W	NovaLife GX 1 Closed Convex	23
932-37	3930W	NovaLife GX 1 Closed Convex	23
932-45	3930W	NovaLife GX 1 Closed Convex	23
941-20	3932Y	NovaLife GX 1 Open Convex	28
941-25	3932Y	NovaLife GX 1 Open Convex	28
941-30	3932Y	NovaLife GX 1 Open Convex	28
941-35	3932Y	NovaLife GX 1 Open Convex	28
941-37	3932Y	NovaLife GX 1 Open Convex	28

Index

NovaLife GX & GX+

Code	SAS Code	Description	Page
941-40	3932Y	NovaLife GX 1 Open Convex	28
941-46	3932Y	NovaLife GX 1 Open Convex	28
941-50	3932Y	NovaLife GX 1 Open Convex	28
941-59	3932Y	NovaLife GX 1 Open Convex	28
942-30	3932Y	NovaLife GX 1 Open Convex	28
942-35	3932Y	NovaLife GX 1 Open Convex	28
942-37	3932Y	NovaLife GX 1 Open Convex	28
942-40	3932Y	NovaLife GX 1 Open Convex	28
942-45	3932Y	NovaLife GX 1 Open Convex	28
942-50	3932Y	NovaLife GX 1 Open Convex	28
942-59	3932Y	NovaLife GX 1 Open Convex	28
952-59	80048C	NovaLife GX+ 1 Open Convex	29
962-10	80069E	NovaLife GX+ 1 Open Flat	25
971-24	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-25	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-30	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-34	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-35	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-40	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-44	80001N	NovaLife GX 1 Closed Soft Convex.....	22
971-54	80001N	NovaLife GX 1 Closed Soft Convex.....	22
981-24	80000M	NovaLife GX 1 Open Soft Convex	26
981-25	80000M	NovaLife GX 1 Open Soft Convex	26
981-30	80000M	NovaLife GX 1 Open Soft Convex	26
981-34	80000M	NovaLife GX 1 Open Soft Convex	26
981-35	80000M	NovaLife GX 1 Open Soft Convex	26
981-40	80000M	NovaLife GX 1 Open Soft Convex	27
981-44	80000M	NovaLife GX 1 Open Soft Convex	27
981-45	80000M	NovaLife GX 1 Open Soft Convex	27
981-54	80000M	NovaLife GX 1 Open Soft Convex	27
982-40	80000M	NovaLife GX 1 Open Soft Convex	27
982-44	80000M	NovaLife GX 1 Open Soft Convex	27
982-45	80000M	NovaLife GX 1 Open Soft Convex	27
984-24	80000M	NovaLife GX 1 Open Soft Convex	27
984-25	80000M	NovaLife GX 1 Open Soft Convex	27
984-30	80000M	NovaLife GX 1 Open Soft Convex	27
984-34	80000M	NovaLife GX 1 Open Soft Convex	27
984-35	80000M	NovaLife GX 1 Open Soft Convex	27
983-24	80000M	NovaLife GX 1 Open Soft Convex	27
983-25	80000M	NovaLife GX 1 Open Soft Convex	27

Index

NovaLife GX & GX+

Code	SAS Code	Description	Page
983-30	80000M	NovaLife GX 1 Open Soft Convex.....	27
983-34	80000M	NovaLife GX 1 Open Soft Convex.....	27
983-35	80000M	NovaLife GX 1 Open Soft Convex.....	27
1301-36	9800T	NovaLife 2 Closed Pouch	30
1301-43	9800T	NovaLife 2 Closed Pouch	30
1301-55	9800T	NovaLife 2 Closed Pouch	30
1301-70	9800T	NovaLife 2 Closed Pouch	30
1302-43	9800T	NovaLife 2 Closed Pouch	30
1302-55	9800T	NovaLife 2 Closed Pouch	30
1302-70	9800T	NovaLife 2 Closed Pouch	30
1303-36	9800T	NovaLife 2 Closed Pouch	30
1303-43	9800T	NovaLife 2 Closed Pouch	30
1303-55	9800T	NovaLife 2 Closed Pouch	30
1315-36	3942L	NovaLife 2 Open Pouch.....	31
1315-43	3942L	NovaLife 2 Open Pouch.....	31
1315-55	3942L	NovaLife 2 Open Pouch.....	31
1315-70	3942L	NovaLife 2 Open Pouch.....	31
1316-36	3942L	NovaLife 2 Open Pouch.....	31
1316-43	3942L	NovaLife 2 Open Pouch.....	31
1316-55	3942L	NovaLife 2 Open Pouch.....	31
1316-70	3942L	NovaLife 2 Open Pouch.....	31
1319-36	3942L	NovaLife 2 Open Pouch.....	31
1319-43	3942L	NovaLife 2 Open Pouch.....	31
1319-55	3942L	NovaLife 2 Open Pouch.....	31
1320-36	3942L	NovaLife 2 Open Pouch.....	31
1321-43	3942L	NovaLife 2 Open Pouch.....	31
1321-55	3942L	NovaLife 2 Open Pouch.....	31
1321-70	3942L	NovaLife 2 Open Pouch.....	31
1322-43	3942L	NovaLife 2 Open Pouch.....	31
1322-55	3942L	NovaLife 2 Open Pouch.....	31
1322-70	3942L	NovaLife 2 Open Pouch.....	31
1736-10	5671L	NovaLife GX 2 Flat Wafer	32
1736-25	5671L	NovaLife GX 2 Flat Wafer	32
1743-10	5671L	NovaLife GX 2 Flat Wafer	32
1743-25	5671L	NovaLife GX 2 Flat Wafer	32
1743-30	5671L	NovaLife GX 2 Flat Wafer	32
1755-10	5671L	NovaLife GX 2 Flat Wafer	32
1755-35	5671L	NovaLife GX 2 Flat Wafer	32
1755-40	5671L	NovaLife GX 2 Flat Wafer	32

Index

NovaLife GX & GX+

Code	SAS Code	Description	Page
1770-10	5671L	NovaLife GX 2 Flat Wafer.....	32
1843-10	80049D	NovaLife GX+ 2 Flat Wafer.....	33
1843-25	80049D	NovaLife GX+ 2 Flat Wafer.....	33
1855-10	80049D	NovaLife GX+ 2 Flat Wafer.....	33
1870-10	80049D	NovaLife GX+ 2 Flat Wafer.....	33
1870-45	80049D	NovaLife GX+ 2 Flat Wafer.....	33
1936-10	3941K	NovaLife GX 2 Convex Wafer.....	34
1936-20	3941K	NovaLife GX 2 Convex Wafer.....	34
1936-23	3941K	NovaLife GX 2 Convex Wafer.....	34
1943-10	3941K	NovaLife GX 2 Convex Wafer.....	34
1943-25	3941K	NovaLife GX 2 Convex Wafer.....	34
1943-30	3941K	NovaLife GX 2 Convex Wafer.....	34
1955-10	3941K	NovaLife GX 2 Convex Wafer.....	34
1955-35	3941K	NovaLife GX 2 Convex Wafer.....	34
1955-40	3941K	NovaLife GX 2 Convex Wafer.....	34
1970-10	3941K	NovaLife GX 2 Convex Wafer.....	34
1970-40	3941K	NovaLife GX 2 Convex Wafer.....	34
1970-45	3941K	NovaLife GX 2 Convex Wafer.....	34
2036-20	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2043-10	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2043-25	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2043-30	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2055-10	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2055-35	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2070-10	80050E	NovaLife GX+ 2 Convex Wafer.....	35
2070-45	80050E	NovaLife GX+ 2 Convex Wafer.....	35
9010-10	80027Y	NovaLife GX 1 Closed Flat.....	21
9010-25	80027Y	NovaLife GX 1 Closed Flat.....	21
9010-30	80027Y	NovaLife GX 1 Closed Flat.....	21
9010-35	80027Y	NovaLife GX 1 Closed Flat.....	21
9010-40	80027Y	NovaLife GX 1 Closed Flat.....	21
9010-45	80027Y	NovaLife GX 1 Closed Flat.....	21
9010-50	80027Y	NovaLife GX 1 Closed Flat.....	21
9050-10	80027Y	NovaLife GX 1 Closed Flat.....	21
9050-25	80027Y	NovaLife GX 1 Closed Flat.....	21
9050-30	80027Y	NovaLife GX 1 Closed Flat.....	21
9050-35	80027Y	NovaLife GX 1 Closed Flat.....	21
9310-25	80017K	NovaLife GX 1 Closed Convex.....	23
9310-30	80017K	NovaLife GX 1 Closed Convex.....	23
9310-35	80017K	NovaLife GX 1 Closed Convex.....	23

Index

NovaLife GX & GX+

Code	SAS Code	Description	Page
9310-50	80017K	NovaLife GX 1 Closed Convex.....	23
9510-20	80048C	NovaLife GX+ 1 Open Convex.....	29
9510-25	80048C	NovaLife GX+ 1 Open Convex.....	29
9510-30	80048C	NovaLife GX+ 1 Open Convex.....	29
9510-35	80048C	NovaLife GX+ 1 Open Convex.....	29
9510-50	80048C	NovaLife GX+ 1 Open Convex.....	29
9510-59	80048C	NovaLife GX+ 1 Open Convex.....	29
9610-10	80069E	NovaLife GX+ 1 Open Flat.....	25
9610-25	80069E	NovaLife GX+ 1 Open Flat.....	25
9610-30	80069E	NovaLife GX+ 1 Open Flat.....	25
9610-35	80069E	NovaLife GX+ 1 Open Flat.....	25
9630-10	80047B	NovaLife GX+ 1 Open Flat.....	25
9630-25	80047B	NovaLife GX+ 1 Open Flat.....	25
9630-30	80047B	NovaLife GX+ 1 Open Flat.....	25

Index

Nova

Code	SAS Code	Description	Page
400-50	N/A	Urostomy Adaptor	50-52
801-20	3802D	Nova 1 Closed Flat	39
801-25	3803E	Nova 1 Closed Flat	39
801-30	3803E	Nova 1 Closed Flat	39
801-35	3803E	Nova 1 Closed Flat	39
801-40	3803E	Nova 1 Closed Flat	39
801-45	3803E	Nova 1 Closed Flat	39
801-50	3803E	Nova 1 Closed Flat	39
802-20	3802D	Nova 1 Closed Flat	39
802-25	3803E	Nova 1 Closed Flat	39
802-30	3803E	Nova 1 Closed Flat	39
802-35	3803E	Nova 1 Closed Flat	39
802-40	3803E	Nova 1 Closed Flat	39
802-45	3803E	Nova 1 Closed Flat	39
802-50	3803E	Nova 1 Closed Flat	39
803-20	3801C	Nova 1 Closed Flat	39
803-25	3801C	Nova 1 Closed Flat	39
803-30	3801C	Nova 1 Closed Flat	39
803-35	3801C	Nova 1 Closed Flat	39
803-40	3801C	Nova 1 Closed Flat	39
805-20	9917Y	Nova 1 Closed Flat	39
806-20	9917Y	Nova 1 Closed Flat	39
815-15	9996D	Nova 1 FoldUp Flat.....	43
816-15	9996D	Nova 1 FoldUp Flat.....	43
818-00	9916X	Nova 1 Infant / Paediatric	47
818-10	9916X	Nova 1 Infant / Paediatric	47
823-15	3873W	Nova 1 FoldUp Flat.....	43
823-20	3873W	Nova 1 FoldUp Flat.....	43
823-25	3873W	Nova 1 FoldUp Flat.....	43
823-30	3873W	Nova 1 FoldUp Flat.....	43
823-35	3873W	Nova 1 FoldUp Flat.....	43
823-40	3873W	Nova 1 FoldUp Flat.....	43
823-45	3873W	Nova 1 FoldUp Flat.....	43
824-15	3873W	Nova 1 FoldUp Flat.....	43
825-15	9947M	Nova 1 FoldUp Flat.....	43
826-15	9947M	Nova 1 FoldUp Flat.....	43
829-20	3871R	Nova 1 Closed Flat	39
829-25	3871R	Nova 1 Closed Flat	39
829-30	3871R	Nova 1 Closed Flat	39
829-40	3871R	Nova 1 Closed Flat	39

Index

Nova

Code	SAS Code	Description	Page
829-50	3871R	Nova 1 Closed Flat	39
831-25	9945K	Nova 1 Closed Convex.....	41
831-30	9945K	Nova 1 Closed Convex.....	41
831-35	9945K	Nova 1 Closed Convex.....	41
831-37	9945K	Nova 1 Closed Convex.....	41
831-40	9945K	Nova 1 Closed Convex.....	41
831-45	9945K	Nova 1 Closed Convex.....	41
832-37	9945K	Nova 1 Closed Convex.....	41
832-46	9945K	Nova 1 Closed Convex.....	41
841-20	9948N	Nova 1 FoldUp Convex.....	46
841-24	9948N	Nova 1 FoldUp Convex.....	46
841-25	9948N	Nova 1 FoldUp Convex.....	46
841-30	9948N	Nova 1 FoldUp Convex.....	46
841-35	9948N	Nova 1 FoldUp Convex.....	46
841-37	9948N	Nova 1 FoldUp Convex.....	46
842-24	9948N	Nova 1 FoldUp Convex.....	46
842-25	9948N	Nova 1 FoldUp Convex.....	46
842-30	9948N	Nova 1 FoldUp Convex.....	46
842-35	9948N	Nova 1 FoldUp Convex.....	46
842-37	9948N	Nova 1 FoldUp Convex.....	46
842-46	9948N	Nova 1 FoldUp Convex.....	46
861-15	9948N	Nova 1 FoldUp X3	44
861-25	9948N	Nova 1 FoldUp X3	44
861-30	9948N	Nova 1 FoldUp X3	44
861-35	9948N	Nova 1 FoldUp X3	44
862-15	9948N	Nova 1 FoldUp X3	44
871-25	9945K	Nova 1 Closed Soft Mouldable	40
871-30	9945K	Nova 1 Closed Soft Mouldable	40
871-35	9945K	Nova 1 Closed Soft Mouldable	40
881-25	9948N	Nova 1 FoldUp Soft Mouldable.....	45
881-30	9948N	Nova 1 FoldUp Soft Mouldable.....	45
881-35	9948N	Nova 1 FoldUp Soft Mouldable.....	45
882-25	9948N	Nova 1 FoldUp Soft Mouldable.....	45
882-30	9948N	Nova 1 FoldUp Soft Mouldable.....	45
882-35	9948N	Nova 1 FoldUp Soft Mouldable.....	45
891-12	9918B	Nova 1 Urostomy X3.....	50
891-20	9918B	Nova 1 Urostomy X3.....	50
891-25	9918B	Nova 1 Urostomy X3.....	50
891-30	9918B	Nova 1 Urostomy X3.....	50
892-12	9918B	Nova 1 Urostomy X3.....	50

Index

Nova

Code	SAS Code	Description	Page
894-20	3837Y	Nova 1 Urostomy Convex.....	51
894-24	3837Y	Nova 1 Urostomy Convex.....	51
894-25	3837Y	Nova 1 Urostomy Convex.....	51
894-30	3837Y	Nova 1 Urostomy Convex.....	51
894-37	3837Y	Nova 1 Urostomy Convex.....	51
894-46	3837Y	Nova 1 Urostomy Convex.....	51
1136-20	3865K	Nova 2 Flat Wafer.....	53
1136-25	3865K	Nova 2 Flat Wafer.....	53
1143-15	3865K	Nova 2 Flat Wafer.....	53
1143-25	3865K	Nova 2 Flat Wafer.....	53
1143-30	3865K	Nova 2 Flat Wafer.....	53
1155-15	3865K	Nova 2 Flat Wafer.....	53
1155-35	3865K	Nova 2 Flat Wafer.....	53
1155-40	3865K	Nova 2 Flat Wafer.....	53
1170-15	3865K	Nova 2 Flat Wafer.....	53
1170-45	3865K	Nova 2 Flat Wafer.....	53
1170-50	3865K	Nova 2 Flat Wafer.....	53
1201-36	3868N	Nova 2 Closed Pouch.....	42
1201-43	3868N	Nova 2 Closed Pouch.....	42
1201-55	3868N	Nova 2 Closed Pouch.....	42
1201-70	3868N	Nova 2 Closed Pouch.....	42
1202-43	3868N	Nova 2 Closed Pouch.....	42
1202-55	3868N	Nova 2 Closed Pouch.....	42
1203-43	9946L	Nova 2 Closed Pouch.....	42
1203-55	9946L	Nova 2 Closed Pouch.....	42
1207-43	80002P	Nova 2 High Output Pouch.....	49
1207-55	80002P	Nova 2 High Output Pouch.....	49
1207-70	80002P	Nova 2 High Output Pouch.....	49
1208-43	80002P	Nova 2 High Output Pouch.....	49
1208-55	80002P	Nova 2 High Output Pouch.....	49
1208-70	80002P	Nova 2 High Output Pouch.....	49
1215-36	3874X	Nova 2 FoldUp Pouch.....	48
1215-43	3874X	Nova 2 FoldUp Pouch.....	48
1215-55	3874X	Nova 2 FoldUp Pouch.....	48
1215-70	3874X	Nova 2 FoldUp Pouch.....	48
1216-36	3874X	Nova 2 FoldUp Pouch.....	48
1216-43	3874X	Nova 2 FoldUp Pouch.....	48
1216-55	3874X	Nova 2 FoldUp Pouch.....	48
1216-70	3874X	Nova 2 FoldUp Pouch.....	48
1217-36	9919C	Nova 2 Urostomy Pouch.....	52

Index

Nova

Code	SAS Code	Description	Page
1217-43	9919C	Nova 2 Urostomy Pouch.....	52
1217-55	9919C	Nova 2 Urostomy Pouch.....	52
1218-36	9919C	Nova 2 Urostomy Pouch.....	52
1218-43	9919C	Nova 2 Urostomy Pouch.....	52
1218-55	9919C	Nova 2 Urostomy Pouch.....	52
1219-36	9949P	Nova 2 FoldUp Pouch.....	48
1219-43	9949P	Nova 2 FoldUp Pouch.....	48
1219-55	9949P	Nova 2 FoldUp Pouch.....	48
1221-55	9852M	Nova 2 FoldUp Pouch.....	48
1221-70	9852M	Nova 2 FoldUp Pouch.....	48
1222-55	9852M	Nova 2 FoldUp Pouch.....	48
1222-70	9852M	Nova 2 FoldUp Pouch.....	48
1343-15	3867M	Nova 2 X3 Wafer.....	54
1343-25	3867M	Nova 2 X3 Wafer.....	54
1355-15	3867M	Nova 2 X3 Wafer.....	54
1355-32	3867M	Nova 2 X3 Wafer.....	54
1355-35	3867M	Nova 2 X3 Wafer.....	54
1443-18	3867M	Nova 2 Soft Mouldable Wafer.....	55
1443-21	3867M	Nova 2 Soft Mouldable Wafer.....	55
1443-25	3867M	Nova 2 Soft Mouldable Wafer.....	55
1443-28	3867M	Nova 2 Soft Mouldable Wafer.....	55
1455-32	3867M	Nova 2 Soft Mouldable Wafer.....	55
1455-35	3867M	Nova 2 Soft Mouldable Wafer.....	55
1455-38	3867M	Nova 2 Soft Mouldable Wafer.....	55
1536-15	3866L	Nova 2 Convex Wafer.....	56
1536-20	3866L	Nova 2 Convex Wafer.....	56
1536-23	3866L	Nova 2 Convex Wafer.....	56
1543-15	3866L	Nova 2 Convex Wafer.....	56
1543-25	3866L	Nova 2 Convex Wafer.....	56
1543-30	3866L	Nova 2 Convex Wafer.....	56
1555-15	3866L	Nova 2 Convex Wafer.....	56
1555-35	3866L	Nova 2 Convex Wafer.....	56
1555-40	3866L	Nova 2 Convex Wafer.....	56
070-20	3879E	Dansac Seal.....	59

Index

Accessories

Code	SAS Code	Description	Page
070-30	3879E	Dansac Seal	59
070-40	3879E	Dansac Seal	59
070-50	3879E	Dansac Seal	59
071-30	80203F	Dansac TRE Seal 30 mm	59
071-40	80203F	Dansac TRE Seal 40 mm	59
072-48	80162C	Dansac TRE Seal 18 mm	59
075-30	80138T	Dansac X-tra Strips.....	60
080-00	3798X	Dansac Nodor S	62
083-01	80051F	Dansac EasiSpray Adhesive Remover	60
085-00	3829M	Dansac Skin Creme.....	61
090-02	3890R	Beige Ostomy Belt.....	62
090-03	3890R	Beige Ostomy Belt.....	62
092-00	3859D	Nova Drainable Clamp.....	63
725-20/30	3539G	Dansac GX-tra Seal.....	59
725-30/30	3539G	Dansac GX-tra Seal.....	59
725-40/30	3567R	Dansac GX-tra Seal.....	59
725-50/30	3567R	Dansac GX-tra Seal.....	59
900-20	3968W	Dansac Complete Irrigation Set (2-piece system)	63
900-35	3968W	Dansac Complete Irrigation Set.....	63
950-20	3916D	Irrigation Sleeves.....	63
950-35	3916D	Irrigation Sleeves.....	63
70000-0000	3716N	Dansac Skin Lotion 200 ml.....	61
70005-0001	N/A	Dansac Skin Lotion 50 ml.....	63
71000-0000	3568T	Dansac Skin Lotion Tissue (Individual Wipes).....	61
71000-0010	9981H	Dansac Skin Lotion Tissue Travel Pack.....	61
77550-0	3552Y	Dansac Soft Paste	60
95200-0000	9804B	Water container	63
95205-0000	3914B	Irrigation Cone with tubing.....	63

Association Contact List

The contact details are correct at time of printing but are subject to change.

Australian Capital Territory

ACT & District Stoma Assn Inc
T/F: (02) 6205 1055
E: stoma@actstoma.onmicrosoft.com

New South Wales

NSW Stoma Limited
T: 1300 OSTOMY
OR (02) 9565 4315
E: info@nswstoma.com.au

Ostomy NSW Ltd
T: (02) 9542 1300
E: orders@ostomynsw.org.au

Northern Territory

Cancer Council of the
Northern Territory
T: (08) 8944 1800
E: ostomy@cancernt.org.au

South Australia

Ostomy Assn of SA
T: (08) 8235 2727
E: colosa@colostomysa.org.au

Ileostomy Assn of SA
T: (08) 8234 2678
E: info@ileosa.org.au

Victoria

Bendigo & District
Ostomy Assn
T: (03) 5441 7520

Colostomy Assn of Vic
T: (03) 9650 1666
E: info@colovic.org.au

Geelong Ostomy Assn
T: (03) 5243 3664
E: goinc@geelongostomy.com.au

Ostomy Assn of Melbourne
T: (03) 9888 8523
E: enquiries@oam.org.au

Peninsula Ostomy Assn
T: (03) 9783 6473
E: poainc1@bigpond.com

Victorian Children's
Ostomy Assn
T: (03) 9345 5325
E: edc@rch.org.au

Warrnambool Ostomy Assn
T: (03) 5563 1446
E: warrnamboolostomy@swh.net.au

Tasmania

Ostomy Tasmania
T: (03) 6228 0799
E: admin@ostomytas.com.au

Queensland

Gold Coast Ostomy Assn
T: (07) 5594 7633
E: gcoa@bigpond.com

North QLD Ostomy Assn
T/F: (07) 4775 2303
E: admin@nqostomy.org.au

QLD Ostomy Assn
T: (07) 3848 7178
E: admin@qldostomy.org.au

QLD Stoma Assn
T: (07) 3359 7570
E: admin@qldstoma.asn.au

Toowoomba and South
West Ostomy Assn
T: (07) 4636 9701
E: bob.schull@bigpond.com

Wide Bay Ostomates Assn
T: (07) 4152 4715
E: wboostomy@bigpond.com

Western Australia

WA Ostomy Assn
T: (08) 9272 1833
E: admin@waostomy.org.au
orders@waostomy.org.au

Terms and Conditions of Sale

LIBERTY MEDICAL PTY LTD TERMS AND CONDITIONS OF SALE

Unless otherwise agreed in writing, these terms and conditions (Conditions) govern all contracts for the sale of goods (Goods) by Liberty Medical Pty Limited (ABN 36 072 817 755) (Company) to the customer (Customer). All other terms and conditions, express or implied, are excluded to the fullest extent permitted by law, including any terms and conditions in the Customer's order for Goods (Order).

1. Quotations, Orders and creation of Contract

Any quotation from Company to Customer is submitted on a confidential basis and is an invitation to treat only, not a contractual offer. Company may, in its absolute discretion, reject any Order or accept any Order in whole or in part.

All orders are subject to acceptance of Liberty Medical and Conditions of Sale. Liberty Medical shall not be liable for any failure to perform the contract of supply of any materials or for any delay or non-delivery due to strikes, fires, explosions, flood, riot, lockouts, injunction, interruption of transportation, accidents, inability to obtain supplies at reasonable prices, war or apprehension of war, force majeure, government action, or any other circumstances beyond Liberty Medical's control.

Liberty Medical reserves the rights to request a minimum order value of \$300 per order. All orders less than \$300 will attract a standard freight charge of \$27.50 + GST.

A contract for the sale of Goods (Contract) will be formed, in accordance with the Order and these Conditions, when Company commences the necessary actions on its part to supply the Goods ordered by Customer.

2. Prices

Goods are invoiced in accordance with the relevant quotation or agreement between Company and Customer, or else from Company's price lists current at the date of delivery.

Prices are quoted exclusive of GST and Customer agrees to pay GST in accordance with the tax invoice.

Liberty Medical endeavours to advise the trade of fluctuations in prices, however all prices are subject to change without notice by Liberty Medical.

3. Payment

Unless otherwise agreed in writing, payment for Goods is due within 30 days of the end of the month in which an invoice is issued without deduction or set off. All payments are to be made in cash, by cheque, credit card or direct deposit into the account notified by Liberty Medical. Failure to pay an amount on the due date entitles Company to, without prejudice to any other remedy, suspend delivery of an Order or cancel any Contract. Customer will be liable for all costs and expenses incurred by Company in obtaining (or attempting to obtain) a remedy for any failure by Customer to pay on the due date.

4. Revision of Terms of Payment

If Liberty Medical at any time considers that conditions do not warrant delivery of the goods on the terms of payment stated herein, Liberty Medical reserves the right to change the term of payment and/or require payment as a prerequisite of delivery and will notify the customer accordingly.

5. Handling

Liberty Medical will select the mode and rate for transport of goods. When special handling/order is requested the Customer will be charged the additional premium freight service, the cost of the special service will be added to the invoice.

6. Risk and delivery

Risk in the Goods will pass to Customer upon delivery. Delivery or dispatch dates are given or accepted by Company in good faith, but not guaranteed unless stated to be "guaranteed" in writing. Company will not be liable for any claim for loss or damages if Company is delayed or prevented from delivering goods for any reason.

7. Title

- a) Property in the goods shall remain with the Company and the Company has clear title to the goods until such time as:
 - (i) full payment is made for all amounts owing by the Customer to the Company so that the Customer's total indebtedness to the Company under these Conditions of Sale is discharged: or
 - (ii) the Customer sells the goods to his customers in the ordinary course of their business.

- b) If the Customer fails to pay any amount of the Customer's indebtedness to the Company under these Conditions of Sale when it is due to the Company or, an event of default as specified in subsection (d) hereof occurs, the Company may withhold notice and without prejudice to any of its other rights and remedies, recover and/or resell the goods or any of them and may enter upon the Customer's premises by its servants or agents for that purpose.
- c) Until full payment the Customer shall store the goods in a way that clearly manifests the Company's title.
- d) Notwithstanding the provisions of these Conditions of Sale payment shall become due immediately upon the Customer committing any act of bankruptcy or the Customer (being a company) committing any act which entitles any person to apply to wind up the Customer, or a Receiver or Receiver and Manager of the Customer is appointed.
- e) The Customer acknowledges that until his total indebtedness under these Conditions of Sale to the Company is discharged he holds the good as bailee of the Company and that a fiduciary relationship exists between them.
- f) In the event of sale of goods to his customers, the Customer in his position as fiduciary shall:
 - (i) assign to the Company the benefit of any claim against a customer, and
 - (ii) account fully to the Company for the proceeds of the sale (unless and until the Customer's total indebtedness to the Company under these Conditions of Sale is discharged).
- g) Until ownership passes to Customer, the Customer:
 - (i) must insure the Goods against all usual risks to full replacement value and hold on trust for Company any insurance monies received for the Goods,
 - (ii) must keep the Goods separate and clearly identified as Company's property,
 - (iii) may sell or otherwise dispose of the Goods in the ordinary course of business, provided Customer holds any proceeds on trust for Company.

8. Returned of Goods

Liberty Medical will accept returns of certain merchandise from the Customer under the following conditions:

- a) All returns must be arranged by Liberty Medical Customer Service Department and a RGA number will be provided. Goods sent without approval will not be accepted and will be re-directed back without issue of credit, at Customer's expense.
- b) Returns must be made within 14 days from the invoice date.
- c) Only full boxes of Goods will be accepted for credit.
- d) Goods accepted for credit where there has been no fault on the part of Liberty Medical, will be charged a re-stocking fee of \$25.00 or 25% of the total value of the returned Goods (whichever is the higher).
- e) Any returned Goods must have an adequate shelf life, require no repackaging and are in an 'as-new' saleable condition. Damaged Goods will not be accepted and will be re-directed back without issue of credit, at Customer's expense.
- f) Freight for goods returned due to Customer's error will be at the Customer's expense.

9. Non-Returnable Merchandise

- a) Sterile Products (except as a result of the Company's shipping error), or by arrangement.
- b) Discontinued items.

10. Product Liability

The product liability of Liberty Medical in relation to any occurrence shall never exceed the limit of Liberty Medical's current product liability insurance.

11. Limitation of Liability

Liberty Medical's liability for any claim whatsoever in respect of products sold hereunder shall not in any case exceed the price of the product in respect to which the claim is made and Customer assumes all risks and liability for the products purchased from Liberty Medical and their use upon and after delivery

including their effect upon any person.

Any and all warranties or conditions express, implied or statutory are excluded except those which cannot be lawfully excluded.

12. Performance

Customer acknowledges, (a) it has relied on its own enquiries concerning the performance of Goods and not on any information, services or recommendation given by Company, and (b) that all descriptive specifications, data, dimensions and weights provided by Company or otherwise contained in Company's catalogues, price lists or other advertising material are approximate only and shall not form part of the contractual description of the Goods unless agreed to in writing by Company, in which case such specifications shall be subject to recognised trade tolerances.

13. Claims

Customer must inspect the Goods upon delivery to ensure the Goods have been delivered in full, are in good working order and free from damage and defects. Any claim by Customer of Goods damaged, lost in transit or short delivered must be made in writing to Company within 7 days of delivery (and must quote the invoice number). Customer must give notice of any defect in the Goods within 7 days of learning of the defect. All Goods that are alleged to be defective must be returned by Customer to Company, if requested, by Company's nominated carrier. Customer will bear the freight cost for the return of Goods which are found not to be defective.

14. Disclaimer

To the extent permitted by law, Company (and its officers, directors, employees, agents or successors) will not be liable for any loss, damage or expense sustained or incurred by Customer or any other party in consequence of or resulting directly or indirectly out of the supply of Goods by Company, the use or receipt of Goods, any breach by Company of any Contract, or the negligence of Company. It is understood that Customer will test and use the Goods as they are intended to be used according to the practices standard in the industry and in strict compliance with all applicable laws and regulations.

Without limiting the generality of the above, Company will not, in any circumstances (including any fault or negligence of Company), be liable for any indirect or consequential losses (including loss of goodwill, business or anticipated savings), loss of profits or use, any rectification costs or any third party claims in connection with Goods or a Contract.

15. Confidentiality

Customer acknowledges that the price of Goods set out in any quotation or Contract, any variation of the payment terms set out in clause 3 and any special conditions agreed in writing between Company and Customer are Confidential Information of Company. Customer may not disclose or otherwise make available such Confidential Information other than to their own employees involved in the Contract, except with Company's prior written consent.

16. Governing Law

Each Contract is governed by the laws of Victoria, Australia and both parties submit to the non-exclusive jurisdiction of the courts of that State.

Liberty Medical Pty Ltd
Trading as Hollister and Dansac
Ground Floor 990 Whitehorse Road
Box Hill, Victoria, Australia 3128

Customer Service
Phone: 1800 279 179
Fax: 03 9899 7587
Email: sales@hollister.com.au
sales@dansac.com.au

P01-550-300 08/19 © 2021 Dansac A/S. AUD049. Updated July 2021.

Dansac Australia
PO Box 375
Box Hill, 3128
Victoria, Australia
Freecall: 1800 880 840

