



Stoma
Skin
Soul



A gentle push
in the right
direction



Introducing the NEW NovaLife TRE™ 2 Piece Soft Convex Barrier

A flexible and adaptable barrier that adheres to your patients' body the way you need it to

- Provides gentle pressure where it's needed, for extra protection against leakage
- Soft and flexible barrier to help you feel comfortable, yet fully supported
- 3 levels of skin health protection with TRE technology
Adhesion, Absorption and pH balance

Comfort, flexibility and skin health combined.
Choose healthy skin for your patients today.

dansac 

The NEW NovaLife TRE™ 2 Piece Soft Convex Barrier

provides all the benefits of a 1 Piece TRE Soft Convex pouch, but with added versatility to meet your patients lifestyle needs.

The benefits of 2 piece convexity

Flexibility

Designed to be flexible to allow for natural body movements

The barrier is only 0.7mm thick compared to a traditional 1.1mm hydrocolloid, the barrier tapers to 0.3mm, making the barrier secure and feel lighter to wear

Softness

TRE barrier is designed to be easy to remove and gentle on the skin


The easily compressible soft convex skin barrier provides gentle pressure where it's needed for extra protection against leakage



Convexity

6mm soft convexity helps push the stoma into the pouch and even out peristomal skin

The depth of convexity and gentle pressure can help with stoma protrusion



Versatility

Pouches can be swapped without having to change the skin barrier every time

The floating flange may help allow a secure pouch and barrier coupling without pressing on the abdomen



Dansac TRE technology offers three levels of protection



pH balance



Absorption



Adhesion

Find out how you can help to maintain healthy skin for your patients with TRE:

To order a sample – 1800 880 840 (Aus) or 0800 678 669 (NZ)

To learn more – visit www.dansac.com.au (Aus) or www.dansac.co.nz (NZ)

Code	Cut	Ring Size	Hole Size
3343-15	cut to fit	43 mm	15-25 mm
3343-25	precut	43 mm	25 mm
3355-15	cut to fit	55 mm	15-38 mm
3355-30	precut	55 mm	30 mm
3355-35	precut	55 mm	35 mm
3370-15	cut to fit	70 mm	15-52 mm

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Dansac, the Dansac logo, NovaLife and TRE are trademarks of Dansac A/S. © 2022 Dansac A/S. AUD137. June 2022.

